# Instructions for use

## Neoss Implant System

Implants, Abutments and associated  
dental, surgical, technical and restorative  
components

**Neoss AB**

Arvid Wallgrens backe 20  
413 46 Göteborg

Sweden

T +46 (0)31 88 12 80

W [www.neoss.com](http://www.neoss.com)

MD 75473

2797

10538\_21 2022-04

## Precautions and Warnings

**DESCRIPTION** Neoss manufactures dental implant components from biocompatible titanium and gold alloys. Accessory components are produced in gold alloy, polymers, titanium and stainless steel.

**PROCEDURES** For detailed information on the specific Neoss Implant System product you are using, please consult the manual, Neoss Implant System Guidelines (10501), instructions for use; Neoss Guide Kit (11408), Neoss TiBase Abutments and ScanPost (11424), Neoss CoCr and Burnout Abutments (10931), Neoss Esthetic Healing Abutments with ScanPeg (11926) and translated versions.

**INDICATIONS** Dental implants from Neoss are used as the means of anchorage between crowns, bridges or dentures and the surrounding bone in the upper or lower jaws. Such prostheses may range from replacement of a single tooth to an entire arch of bridgework and cement or screw retained restorations. The Neoss Implants – Neoss Implant System are also intended for immediate placement and function on single tooth and/or multiple tooth applications recognizing sufficient bone stability and appropriate occlusal loading, to restore chewing function. Multiple tooth applications may be splinted with a bar.

**CONTRAINdications** Insufficient number of implants or size of implants or abutments to support biomechanical loads or undesirable positioning of implants can lead to mechanical failures including fatigue fracture of implant fixtures, abutments or abutment screws. Such an example is narrow diameter implants in combination with angulated abutments such as angulated Access Abutments in the posterior region. Implant placement and prosthetic design must accommodate individual patient conditions such as bruxism or unfavorable jaw relationships to reduce the risk of overload or fatigue failure, and treatment is contraindicated if adequate accommodation cannot be accomplished. Please consult appropriate surgical and restorative manuals and textbooks for information on treatment planning and medical evaluation.

Dental implants need to be selected so that they provide an optimal connection and stability with the surrounding bone. Neoss dental implants are contraindicated where bone of insufficient quality or quantity is available to produce adequate stability and support. The design and construction of the abutment and prosthesis by the technician should incorporate appropriate retentive features for the prosthesis and should optimize the angulation between the implant fixtures and prosthesis such that applied loads are directed down the long axis of the implant. Failure to achieve this can lead to excessive bending force and fatigue failure of the implant components.

The Neoss Implant System has specific design characteristics for mating Neoss components such as implants, abutments and prosthetic components. Combining components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory esthetic results.

Pre-operative hard tissue or soft tissue deficits may result in a compromised esthetic outcome or unfavorable implant angulation.

**GENERAL PRECAUTIONS** Surgical and restorative products used to achieve and maintain osseointegration should be utilized by persons trained in this method. Such training is offered at a number of centers. Please contact the manufacturer for information. Pre-operative patient evaluation and close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for success. Neoss implants, abutments and abutment screws must be used solely on one patient. Single use devices should not be reused due to risks of product contamination, patient/user infection and/or failure of the device to perform as intended. In the event of malfunction of the device or changes in its performance of the device, the patient should contact the dentist for assessment. Handling of hazardous material according to established procedures at the hospital/clinic.

**PROCEDURAL PRECAUTIONS** Especially important is proper stress distribution: relining of the denture after implant installation to avoid premature loading; passive adaptation and fitting of the bridge to the implant abutments; adjusting occlusion to the opposing jaw; avoiding excessive transverse loading forces are all important considerations. Prosthesis metal substructures should be made using either titanium or suitable gold and CoCr alloys. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

**ADVERSE EFFECTS** Implant techniques have normal contraindications and risks. These are extensively documented in the dental literature. Incorrect clinical placement or loading resulting in loss of implant anchorage or loss of the prosthesis are possible events after surgery. Lack of bone quantity or quality, infections, poor patient hygiene or cooperation, and general diseases are some potential causes for loss of anchorage and function.

**STERILITY** All dental implant, bone cutting instruments and some healing abutments are supplied sterile with a given expiry date as indicated by the packaging. For re-sterilization of non-sterile or sterile instruments and accessories please refer to IFU Cleaning and Maintenance of Neoss Reusable Products (14077).

## LIFETIME

**Implants, permanent abutments, permanent abutment/prosthetic screws:** The patient can, if not belonging to any groups with contraindications and if device is used as intended and for correct indications, expect the dental implant treatment to restore and maintain aesthetics and/or function for a time period comparable to natural teeth.

**Healing abutments, Cover Screws, Provisional abutments and Provisional Abutment/Prosthetic Screws:** The patient can, if not belonging to any groups with contraindications and if the device is used as intended and for correct indications, expect the devices to stay in place until replaced with abutments and restorations.

**Overdenture Attachments:** With repeated snapping in and out of the denture prosthesis, the overdenture attachments eventually wear out over time, resulting in diminished snap retention and perceived looseness of the overdenture. When this happens they are to be replaced with new attachments.

**Bonecutting instruments:** A maximum of 10 re-use cycles are recommended for the multiple use devices. However, the instruments shall always be checked for damage before re-use. If damage is found then the instrument is to be disposed of.

**Instruments, Scanpost, Trays:** The instruments and accessories are reusable devices that can be reused as long as the integrity and performance of the device are maintained. The devices must be inspected before each use for visible signs of damage, deformation, wear or corrosion. Devices showing any signs of visible damage or loss of functional compatibility (for example fit of instrument with mating components, lifting function etc) shall be discarded.

**STORAGE** Sterilized bags and unused components must be stored in dry environment, at room temperature and out of direct sunlight. Unused components must be stored in their original packaging.

**MRI SAFETY INFORMATION** Neoss dental implants and restorative components have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of these devices in the MR environment is unknown. Scanning a patient who has these devices may result in patient injury.

**CAUTION** Federal (USA) law restricts the sale of this device to or on the order of a licensed physician or dentist.

## GENERAL PACKAGING SYMBOLS

USE BY/EXPIRY DATE	CATALOGUE NUMBER	LOT/BATCH NUMBER	DO NOT RE-USE (Single use only)	TEMPERATURE LIMIT
KEEP AWAY FROM SUNLIGHT	MANUFACTURER	DATE OF MANUFACTURE	STERILIZED USING ETHYLENE OXIDE	STERILE BY IRRADIATION (Contents of inner package sterile)
NON-Sterile	DO NOT USE IF PACKAGE IS DAMAGED	MEDICAL DEVICE	CONSULT INSTRUCTIONS FOR USE (Also available on <a href="http://www.neoss.com/IFU">www.neoss.com/IFU</a> )	CAUTION: Federal (USA) law restricts the sale of this device to or on the order of a licensed physician or dentist <b>Rx only</b>

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